

PDREX[®]

Synthetic Absorbable
POLYDIOXANONE
SUTURES



PDREX[®] SUTURES

PDREX[®] sutures is synthetic absorbable monofilament suture composed of a poly (pdioxanone) and is available dyed

When dyed, only FDA –approved color additives such as D&C Violet No. 2 are used.

PDREX[®] meets all the requirements of the European Pharmacopoeia (EP) for Synthetic Absorbable Monofilament Sutures.

PDREX[®] has been found th be non antigenic and elicits only minimal tissue reactivity during the absorption process. PDREX[®] sutures provide wound support for longer periods up th 50 days and offers far superior tensile strength and outstanding pliability. Its monofilament structure provieds good handling properties and excellent knot security. PDREX[®] retains approximately 50~60% of its original tensile strength at 4 weeks post implantation, Complete absorption in tissues takes around 180 to 210 days.

Needle Point Type

- Taper point ○
- Reverse Cutting ▼
- Conventional cutting ▲
- Spatula ▽
- Blunt ⊙
- Taper Cutting ⊕






















PDREX[®]













PDO (Polydioxanone) - Monofilament

Needle Type	Length & Color	Suture U.S.P (Metric)				
		7-0 (0.5)	6-0 (0.7)	5-0 (1)	4-0 (1.5)	3-0 (2)
6mm ▽ 3/8	30cm Violet	PB706				
	30cm Violet	PB706D				
6.5mm ○ 1/2	70cm Violet	PN765				
7mm ▽ 3/8	45cm Violet	PB707				
8mm ▽ 1/2	45cm Violet	PA708				
	45cm Violet	PA708D				
1/4	20cm Violet		PA608SFD			
	45cm Violet		PA608FD	PA508FD	PA408FD	
○ 3/8	30cm Violet		PM608D			
	30cm Violet		PM608S			
	45cm Violet		PM608			
▽ 1/2	45cm Violet			PF508		
9mm ○ 1/2	45cm Violet		PN609			
▽ 3/8	45cm Violet		PE609			
	45cm Violet			PE509D		
9.3mm ○ 3/8	70cm Violet	PM793	PM693			
10mm ○ 1/2	70cm Violet		PN610			
11mm ▽ 3/8	45cm Violet		PE611			

PDREX[®]




















PDO (Polydioxanone) - Monofilament

















Needle Type	Length & Color	Suture U.S.P (Metric)				
		6-0 (0.7)	5-0 (1)	4-0 (1.5)	3-0 (2)	2-0 (3)
11mm						
 5/8	 30cm Violet	PT611				
3/8	 45cm Violet	PT611FD	PT511FD	PT411FD		
 3/8	 70cm Violet	PM611				
13mm						
 3/8	 45cm Violet		PB513			
str.	 45cm Violet			PS413		
	45cm Violet	PM613S				
 3/8	 60cm Violet			PM413		
	70cm Violet	PM613	PM513			
1/2	45cm Violet		PN513S			
	70cm Violet	PN613	PN513	PN413	PN313	PN213
 3/8	 45cm Violet	PE613	PE513	PE413		
 3/8	 70cm Violet	PT613F				
1/2	 70cm Violet		PT513F			
14mm						
 1/2	45cm Violet		PN514S	PN414S		
	70cm Violet	PN614	PN514	PN414		
15mm						
 3/8	 45cm Violet			PM415		
1/2	 70cm Violet		PN515			
 3/8						

Needle Type	Length & Color	Suture U.S.P (Metric)				
		6-0 (0.7)	5-0 (1)	4-0 (1.5)	3-0 (2)	2-0 (3)
16mm						
⊙ 3/8 	75cm Violet			PM416		
1/2 	70cm Violet				PN316	PN216
△ 3/8 	60cm Violet			PEE416S		
	75cm Violet			PEE416		
17mm						
⊙ 3/8 	45cm Violet			PM417		
	90cm Violet		PM517			
1/2 	75cm Violet	PN617	PN517	PN417	PN317	PN217
▽ 3/8 	75cm Violet		PE517	PE417	PE317	
⊕ 1/2 	75cm Violet		PT517	PT417	PT317	PT217
18mm						
▽ 3/8 	75cm Violet			PE418		
⊙ 3/8 	75cm Violet		PM518	PM418	PM318	
	1/2 	75cm Violet	PN618			
19mm						
▽ 3/8 	30cm Violet			PE419S		
	45cm Violet			PE419		
⊙ 1/2 	75cm Violet		PN519	PN419	PN319	
str. 	35cm Violet		PNS519	PN419S		

PDREX[®]


















PDO (Polydioxanone) - Monofilament








Needle Type	Length & Color	Suture U.S.P (Metric)				
		4-0 (1.5)	3-0 (2)	2-0 (3)	0 (3.5)	1 (4)
20mm						
 1/2 	75cm Violet	PN420	PN320			
21mm						
 1/2 	75cm Violet	PN421	PN321			
 1/2 	75cm Violet	PT421	PT321			
 3/8 	75cm Violet		PE321	PE221		
	45cm Violet			PF221		
1/2 	75cm Violet		PF321			
22mm						
 1/2 	75cm Violet	PN422	PN322	PN222	PN022	
 1/2 	75cm Violet	PT422	PT322			
 1/2 	75cm Violet					PF0122
23mm						
 1/2 	75cm Violet		PF323	PF223		
 str. 	30cm Violet				PNS023	

Needle Type	Length & Color	Suture U.S.P (Metric)				
		4-0 (1.5)	3-0 (2)	2-0 (3)	0 (3.5)	1 (4)
24mm						
 3/8	 75cm Violet		PE324	PE224		
 1/2	 75cm Violet	PN424				
25mm						
 3/8	 75cm Violet		PE325	PE225		
1/2	 75cm Violet	PN425	PN325	PN225	PN025	
 1/2	 75cm Violet		PT325	PT225		
26mm						
 1/2	 75cm Violet	PN426	PN326	PN226	PN026	PN0126
	 90cm Violet	PN426L	PN326L	PN226L		
3/8	 75cm Violet		PM326			
5/8	 75cm Violet		PN326F	PN226F		
 1/2	 75cm Violet	PT426	PT326	PT226		

PDREX[®]













PDO (Polydioxanone) - Monofilament






Needle Type	Length & Color	Suture U.S.P (Metric)				
		4-0 (1.5)	3-0 (2)	2-0 (3)	0 (3.5)	1 (4)
26mm						
 3/8 	75cm Violet			PE226		
1/2 	75cm Violet		PF326	PF226	PF026	PF0126
27mm						
 5/8 	75cm Violet			PN227F		
28mm						
 1/2 	75cm Violet		PU328	PU228		
 1/2 	75cm Violet		PN328	PN228	PN028	PN0128
 	75cm Violet			PN228D		
 5/8 	75cm Violet			PN228F	PN028F	
29mm						
 1/2 	75cm Violet		PN329	PN229	PN029	PN0129
 1/2 	75cm Violet			PT229	PT029	PT0129

Needle Type	Length & Color	Suture U.S.P (Metric)				
		4-0 (1.5)	3-0 (2)	2-0 (3)	0 (3.5)	1 (4)
30mm						
⊙ 3/8 	75cm Violet			PE230		
○ 1/2 	75cm Violet	PN430	PN330	PN230	PN030	PN0130
31mm						
⊙ 1/2 	75cm Violet			PN2231S	PN031	PN0131
	150cm Violet			PN231 (LOOP)		
33mm						
⊙ 1/2 	75cm Violet			PU233	PU033	
○ 1/2 	75cm Violet			PN233S	PN033S	
	90cm Violet			PN233	PN033	
35mm						
⊙ 1/2 	75cm Violet		PN335	PN235	PN035	
	90cm Violet					PN0135
36mm						
⊙ 1/2 	75cm Violet	PN436	PN336	PN236	PN036	PN0136
	90cm Violet			PN236L		PN0136L
▽ 1/2 	75cm Violet			PF236	PF036	PF0136
⊕ 1/2 	75cm Violet					PT0136

PDREX[®]











PDO (Polydioxanone) - Monofilament






Needle Type	Length & Color	Suture U.S.P (Metric)				
		3-0 (2)	2-0 (3)	0 (3.5)	1 (4)	2 (5)
37mm						
 3/8 	75cm Violet	PN337	PE237			
 1/2 	75cm Violet		PT237	PT037		
 1/2 	75cm Violet	PN337	PN237	PN037	PN0137	
40mm						
 1/2 	75cm Violet 90cm Violet			PF040	PF0140S PF0140	PF0240
 1/2 	150cm Violet			PT040 (LOOP)		
 1/2 	75cm Violet 90cm Violet		PN240S PN240	PN040S PN040	PN0140S	PN0240

Needle Type	Length & Color	Suture U.S.P (Metric)				
		3-0 (2)	2-0 (3)	0 (3.5)	1 (4)	2 (5)
40mm						
⊙ 5/8 	150cm Violet			PN040L (LOOP)	PN0140 (LOOP)	
41mm						
▽ 1/2 	75cm Violet			PF041	PF0141	
⊕ 1/2 	75cm Violet			PT041	PT0141	
⊙ 1/2 	75cm Violet 90cm Violet			PU041S PU041	PU0141S PU0141	
⊙ 1/2 	75cm Violet 90cm Violet		PN241S PN241	PN041S PN041	PN0141S PN0141	

PDREX[®]









PDO (Polydioxanone) - Monofilament





Needle Type	Length & Color	Suture U.S.P (Metric)				
		3-0 (2)	2-0 (3)	0 (3.5)	1 (4)	2 (5)
45mm						
 3/8 	75cm Violet		PE245	PE045		
 1/2 	90cm Violet					PT0245
 1/2 	90cm Violet				PU0145	
 1/2 	75cm Violet					PN0245S
	90cm Violet				PN0145	PN0245
46mm						
 1/2 	75cm Violet				PT0146	PT0246

Needle Type	Length & Color	Suture U.S.P (Metric)				
		3-0 (2)	2-0 (3)	0 (3.5)	1 (4)	2 (5)
46mm						
⊙ 1/2 	75cm Violet			PN046	PN0146	PN0246
48mm						
▽ 1/2 	75cm Violet					PF0248S
	90cm Violet			PF048	PF0148	PF0248
⊕ 1/2 	150cm Violet			PT048 (LOOP)		
⊙ 1/2 	150cm Violet			PU048 (LOOP)	PU0148 (LOOP)	
	75cm Violet					PN0248
⊙ 1/2 	90cm Violet			PN048S	PN0148S	
	150cm Violet			PN048 (LOOP)	PN0148 (LOOP)	

PDREX[®]

PDO (Polydioxanone) - Monofilament

Needle Type	Length & Color	Suture U.S.P (Metric)				
		3-0 (2)	2-0 (3)	0 (3.5)	1 (4)	2 (5)
49mm						
 1/2 	90cm Violet					PF0249
 1/2 	75cm Violet			PN049	PN0149S	PN0249
	90cm Violet				PN0149	
63mm						
 3/8 	90cm Violet			PU063	PU0163	
65mm						
 1/2 	75cm Violet					PN0265
	150cm Violet			PN065 (LOOP)	PN0165 (LOOP)	

Needle Type	Length & Color	Suture U.S.P (Metric)				
		3-0 (2)	2-0 (3)	0 (3.5)	1 (4)	2 (5)
65mm  str. 	75cm Violet			PNS065		
70mm  1/2 	75cm Violet			PN070		





46, Nakdong-daero 1302beon-gil, Sasang-gu, Busan, Korea
TEL : 82-51-305-8016 FAX : 82-51-305-8021
E-mail : sales@sm-eng.net http://www.sm-eng.net