

A bright white starburst graphic with multiple sharp rays emanating from a central point, positioned above the 'Rex Point' text.

Rex Point®

A large, thick white arc that starts from the left side of the page and curves upwards and to the right, framing the central text.

STAINLESS STEEL
SUTURE NEEDLES





C O N T E N T S



04 | TYPE OF NEEDLE

05 | TYPE OF CIRCLE · TYPE OF EYE

06 | 외과용 각침 3/8 CIRCLE(GENERAL CLOSURE) · 외과용 각침 1/2 CIRCLE(GENERAL CLOSURE)

07 | 외과용 각침 STRAIGHT(GENERAL CLOSURE) · 내복용 환침 2 1/2 CIRCLE(ABDOMINAL)

08 | 내장용 환침 1/2 CIRCLE(INTESTINAL) · 내장용 환침 STRAIGHT(INTESTIANL)

09 | 부인과용 환침 1/2 CIRCLE(OB-GYN) · 부인과용 각침 1/2 CIRCLE(VASCULAR)

10 | 안과용 각침 3/8 CIRCLE(OPHTHALMIC) · 혈관용 환침 1/2 CIRCLE(VASCULAR)

11 | 피부용 각침 3/8 CIRCLE(SKIN) · 각막용 각침 3/8 CIRCLE(CORNEAL)

TYPE OF NEEDLE




TAPER POINT NEEDLE

THIS TYPE OF NEEDLE PENETRATES, PASSES THROUGH TISSUES BY STRECHING WITHOUT TISSUE CUTTING, NORMALLY THIS NEEDLE MINIMIZES POTENTIAL TEARING OF FASCIA AND IS USED FOR EASILY PENETRATED TISSUES AS LIKE ABDOMINAL VISCERA, PERITONEUM, DURA, SUBCUTANEOUS LAYERS.



CONVENTIONAL CUTTING NEEDLE

THIS TYPE OF NEEDLE HAS 3 CUTTING EDGES AND DESIGNED FOR TISSUE THAT IS TOUGH TO PENETRATE AS LIKE SKIN, TENDON SHEATHS, ORAL MUCOSA.



REVERSE CUTTING NEEDLE




SAME WITH ABOVE CONVENTIONAL TYPE AND THESE NEEDLES ARE ALSO HELPFUL IN COSMETIC AND OPHTHALMIC SURGERY, CAUSING IT'S MINIMAL TRAUMA.



SPATULA CUTTING NEEDLE


THIS TYPE OF NEEDLE DESIGNED TO REDUCED TISSUE INJURY AND ALLOW MAXIMUM EASE OF PENETRATION AND CONTROL AS IT PASS BETWEEN AND THROUGH TISSUE LAYERS.

TYPE OF CIRCLE



TYPE OF EYE

REGURAL EYE



SPRING EYE












외과용 각침 3/8 CIRCLE (GENERAL CLOSURE)

P102 : REGULAR EYE





F102 : SPRING EYE



외과용 각침 1/2 CIRCLE (GENERAL CLOSURE)

P106 : REGULAR EYE

F106 : SPRING EYE






외과용 각침 STRAIGHT (GENERAL CLOSURE)

P162 : REGULAR EYE


F162 : SPRING EYE



내복용 환침 1/2 CIRCLE (ABDOMINAL)

P114 : REGULAR EYE

F114 : SPRING EYE






내장용 환침 1/2 CIRCLE (INTESTINAL)

P136 : REGULAR EYE


F136 : SPRING EYE



내장용 환침 1/2 STRAIGHT (INTESTINAL)

P134 : REGULAR EYE

F134 : SPRING EYE











부인과용 환침 1/2 CIRCLE (OB-GYN)

P144 : REGULAR EYE







F144 : SPRING EYE



부인과용 각침 1/2 CIRCLE (OB-GYN)

P142 : REGULAR EYE

F142 : SPRING EYE





혈관용 환침 1/2 CIRCLE
(VASCULAR)

P24 : REGULAR EYE


F24 : SPRING EYE



혈관용 환침 3/8 CIRCLE
(VASCULAR)

P22 : REGULAR EYE

F22 : SPRING EYE





안과용 각침 3/8 CIRCLE (OPHTHALMIC)

P124 : REGULAR EYE

F124 : SPRING EYE

00 11

0 13

1 15

2 18

3 21



각막용 각침 3/8 CIRCLE (CORNEAL)

P128 : REGULAR EYE

F128 : SPRING EYE

1 5

2 6

3 7

4 8


5 9



피부용 각침 3/8 CIRCLE
(SKIN)

P104 : REGULAR EYE









F104 : SPRING EYE



외과용 각침 선단곡형
(GENERAL CLOSURE)

P108 : REGULAR EYE

F108 : SPRING EYE





46, Nakdong-daero 1302beon-gil, Sasang-gu, Busan, Korea
TEL : 82-51-305-8016 FAX : 82-51-305-8021
E-mail : sales@sm-eng.net http://www.sm-eng.net